

Materials List

Zest™

Component	Material
Mask Base	Polycarbonate (PC)
Elbow	Polycarbonate (PC)
Swivel	Polycarbonate (PC)
Silicone Seal	Silicone (LSR) / Polycarbonate (PC)
Forehead Cushion	Silicone (LSR)
Foam Cushion	Polyurethane foam
Glider™ Strap	Nylon (PA)
Stretchgear™ Headgear	Breathe-o-prene™ (Lycra™ / Polyurethane foam laminate)
Pressure Ports Connector (Euro Only)	Polycarbonate (PC) / Santoprene™ (TPE)

Opus™ 360

Component	Material
Mask Base	Polycarbonate (PC)
Ball Elbow	Polycarbonate (PC)
Swivel Snap 17mm	Polycarbonate (PC)
Swivel Snap Male	Polycarbonate (PC)
Flexitube	Polyethylene (PE)
Silicone Seal (S, M, L)	Silicone (LSR)
Rigid Headgear	Nylon (PA)
Loop Tab	Velcro™
Headgear	Breathe-o-prene™ (Lycra™ / Polyurethane foam laminate)
Mask Frame Pad	Softlon 3000 (Polyethylene (PE) foam)

Forma™

Component	Material
Mask Base	Polycarbonate (PC)
Non-Rebreathing Valve	Polycarbonate (PC) / Silicone (LSR)
Elbow	Polycarbonate (PC)
Glider™ Strap	Nylon (PA)
Glider™ Clip	Acetal (POM)
Forehead Cushion	Silicone (LSR)
Silicone Seal (S, M, L, XL)	Silicone (LSR)
Foam Cushion (S, M, L, XL)	Polyurethane foam
Headgear	Breathe-o-prene™ (Lycra™ / Polyurethane foam laminate)
Pressure Ports Connector (Euro Only)	Polycarbonate (PC) / Santoprene™ (TPE)

Materials List

FlexiFit™ 405

Component	Material
Mask Base	Polycarbonate (PC)
Elbow	Polycarbonate (PC)
Exhaust Cover	Polycarbonate (PC)
Exhaust Diffuser	Polyester fiber
Swivel	Polycarbonate (PC)
Silicone Seal (S / L)	Silicone (LSR)
Forehead Cushion	Silicone (LSR)
Foam Cushion	Polyurethane foam
Glider™ Strap	Nylon (PA)
Headgear	Breathe-o-prene™ (Lycra™ / Polyurethane foam laminate)
Pressure Ports Connector (Euro Only)	Polycarbonate (PC) / Santoprene™ (TPE)

FlexiFit™ 406

Component	Material
Mask Base	Polycarbonate (PC)
Elbow	Polycarbonate (PC)
Swivel	Polycarbonate (PC)
Silicone Seal	Silicone (LSR)
Forehead Cushion	Silicone (LSR)
Foam Cushion	Polyurethane foam
Glider™ Strap	Nylon (PA)
Stretchgear™ Headgear	Breathe-o-prene™ (Lycra™ / Polyurethane foam laminate)
Pressure Ports Connector (Euro Only)	Polycarbonate (PC) / Santoprene™ (TPE)

FlexiFit™ 407

Component	Material
Mask Base	Polycarbonate (PC)
Elbow	Polycarbonate (PC)
Swivel	Polycarbonate (PC)
Silicone Seal	Silicone (LSR)
Forehead Cushion	Silicone (LSR)
Foam Cushion	Polyurethane foam
Glider™ Strap	Nylon (PA)
Stretchgear™ Headgear	Breathe-o-prene™ (Lycra™ / Polyurethane foam laminate)
Pressure Ports Connector (Euro Only)	Polycarbonate (PC) / Santoprene™ (TPE)

Materials List

FlexiFit™ 431

Component	Material
Mask Base	Polycarbonate (PC)
Non-Rebreathing Valve	Polycarbonate (PC) / Silicone (LSR)
Elbow	Polycarbonate (PC)
Glider™ Strap	Nylon (PA)
Glider™ Clip	Acetal (POM)
Forehead Cushion	Silicone (LSR)
Silicone Seal (S, M, L)	Silicone (LSR)
Headgear	Breathe-o-prene™ (Lycra™ / Polyurethane foam laminate)
Pressure Ports Connector (Euro Only)	Polycarbonate (PC) / Santoprene™ (TPE)

FlexiFit™ 432

Component	Material
Mask Base	Polycarbonate (PC)
Non-Rebreathing Valve	Polycarbonate (PC) / Silicone (LSR)
Elbow	Polycarbonate (PC)
Glider™ Strap	Nylon (PA)
Glider™ Clip	Acetal (POM)
Forehead Cushion	Silicone (LSR)
Silicone Seal (S, M, L)	Silicone (LSR)
Foam Cushion (S, M, L, XL)	Polyurethane foam
Headgear	Breathe-o-prene™ (Lycra™ / Polyurethane foam laminate)
Pressure Ports Connector (Euro Only)	Polycarbonate (PC) / Santoprene™ (TPE)

Materials List

Oracle™

Component	Material
Mask Base	Polycarbonate (PC) / Silicone (LSR)
Non-Rebreathing Valve	Polycarbonate (PC) / Silicone (LSR)
Elbow	Polycarbonate (PC)
Exhaust Cover	Polycarbonate (PC)
Exhaust Diffuser	Polyester fiber
Male/Female Swivel	Polycarbonate (PC)
Flexitube	Polyethylene (LLDPE)
Slider	Polyetherimide (PEI)
Soft Seal (M / L)	Silicone (LSR)
Nasal Plugs	Silicone (LSR)
Headgear	Breathe-o-prene™ (Lycra™ / Polyurethane foam laminate)
Pressure Ports Connector (Euro Only)	Polycarbonate (PC) / Santoprene™ (TPE)

Aclaim®

Component	Material
Mask Base	Polycarbonate (PC)
T-Piece	Polycarbonate (PC)
Exhaust Cover	Polycarbonate (PC)
Exhaust Diffuser	Polyester fiber
Swivel	Polycarbonate (PC)
Flexitube	Polyethylene (LLDPE)
Flexitube Cuff	Polypropylene (PP)
Silicone Seal (S / L)	Silicone (LSR)
Forehead Cushion	Silicone (LSR)
Foam Cushion	Polyurethane foam
Glider™ Strap	Nylon (PA)
Headgear	Breathe-o-prene™ (Lycra™ / Polyurethane foam laminate)
Pressure Ports Connector (Euro Only)	Polycarbonate (PC) / Santoprene™ (TPE)