

1 inch



- Full Face Delivery
- Direct Nasal Delivery
- Nasal Delivery
- Oral Delivery




**Infinity<sup>481</sup>**  
DIRECT NASAL MASK



**FlexiFit<sup>407</sup>**  
PREMIUM NASAL MASK



**FlexiFit<sup>406</sup>**  
PETITE NASAL MASK



### Fisher & Paykel Multi-Mask Sizing Gauge

1. Print document - Note do not choose to scale document to paper size when printing. (Square at top of page should be one inch wide/tall)
2. Cut out paper inside of mask outlines
3. Choose correct size according to fit (406/407) or size lines (431/432/481).

**FlexiFit<sup>431/432</sup>**  
FULL FACE MASK

185044349REVA



Presented With Compliments From  
[www.DirectHomeMedical.com](http://www.DirectHomeMedical.com)

**Fisher & Paykel**  
HEALTHCARE

[www.fphcare.com](http://www.fphcare.com)